

# KWIK KLEEN®

## KERRISON

### FEATURES:

- Patented design enables the Kerrison to be fully cleaned, thus eliminating the build up of bone, tissue, and bio-burden that can cause sticking and non-performance.
- Wide grip handles provide complete comfort and control.
- Single lever makes ease of opening and closing.
- Bite size numbers are etched on both sides of Kerrisons.
- Proven design provides a consistent smooth cutting action.
- 180 mm - 200 mm - 230 mm length shaft available - also custom lengths available.
- Top rail does not detach, thus eliminating errors in assembly.
- Manufactured from 100% high grade German stainless steel.

