

SPEEDLOCK[®]

Laparoscopic Instruments

SPEEDLOCK[®] - System

In terms of quality and handling these instruments surpass any other product line in the endoscopic market today. Some highlights include a cleaning port, which allows easy and time saving reprocessing, without taking the instruments apart. Jaw made out of special stainless steel, making them even harder and better resistant against torque, tension and corrosion. The innovative new and unique bayonet coupling, makes assembly secure and so much easier.

Only an ongoing dialogue with users in combination with our multiple years of experience and a continuous development have made this exceptional modular instrument line possible and that's why using SPEEDLOCK[®] means being one step ahead.

- Rotatable with custom color rings for rotators
- Demountable (handle, sheath and insert)
- Flushable
- Unique bayonet coupling
- 45 degree handle-top HF-port
- Cable protection
- Length 26cm - 30cm - 36cm - 45cm "no problem with obese patients"
- Ergonomic handles in different styles for everybody's hand, large and small
- 10 different handle styles available
- 360 degree rotation, no increments to distort positioning
- No hidden spaces to accommodate blood and debris
- Time saving - through easy cleaning and assembly
- Highest durability, even if autoclaved
- Cost effective - SPEEDLOCK[®] design reduces cost per case and inventory management
- Virtually no stock if compared to disposable instruments
- Environment - friendly due to its outstanding durability

Interchangeable color coded rings

